

LUNCH & LEARN LIVE WEBINAR

# Mastering MetaXplore™:

Interpret Microbiome Reports in Under 15 Minutes

Wednesday 4 September | 12PM AEST

With Dr Brad Leech



1

## Meet your host



**Dr Brad Leech**  
Nutritionist and Lead Clinical Educator



All participants have been muted



There is an optional 15 minutes for questions at the end



Add your questions in the chat to have them answered live

**CO-BIOME**  
MICROBIOME

2

## Disclaimers

- The information provided in this webinar is for the use of qualified healthcare professionals.
- The information contained in this webinar is in no way to be taken as prescriptive or to replace a healthcare professional's duty of care and personalised care practices.
- The clinical opinions and patient case studies shared by presenters are solely those of the individual presenters and do not necessarily represent the view of Co-Biome.



CO-BIOME  
by MICROBA

3

## Objectives

1. Learn how to interpret gastrointestinal health and the health of the whole gut microbiome in under 15 minutes using the report interpretation checklist
  - ✓ Identify red flags and when to refer
  - ✓ Assess gut terrain
  - ✓ Identify functional dysbiosis
2. Understand how to save time by avoiding common roadblocks when interpreting a MetaXplore™ report
3. Learn how to confidently apply report insights to your patient's management plan
4. Discover how MetaXplore's Expert Summary is helping you interpret patient reports in just minutes



CO-BIOME  
by MICROBA

4

# When to consider assessing the gut microbiome and gastrointestinal health




- Gastrointestinal disorders, including irritable bowel syndrome (IBS), intestinal permeability, and inflammatory bowel disease (IBD)
- Digestive complaints, including constipation, diarrhoea, bloating and abdominal pain
- Hormonal imbalance
- Metabolic and weight concerns
- Immune system health concerns
- Chronic inflammation



5

MetaXplore™ Range Report Interpretation Checklist



1 ASSESS	Red flags (refer to a medical specialist, if necessary)	Faecal occult blood detected*	<input type="checkbox"/>	
		Calprotectin above 100 µg/g*	<input type="checkbox"/>	
		Lactoferrin above 7.2 µg/g*	<input type="checkbox"/>	
		Pancreatic elastase below 100 µg/mL*	<input type="checkbox"/>	
		Pathogens detected on diagnostic targeted pathogen panel**	<input type="checkbox"/>	
	Gut terrain	Potential pathogens identified in metagenomic species table (search pathogen)	<input type="checkbox"/>	
		Faecal pH*	<input type="checkbox"/>	
		Secretory IgA*	<input type="checkbox"/>	
		Zonulin*	<input type="checkbox"/>	
		Mucin degradation	<input type="checkbox"/>	
Dysbiosis	Oral species	<input type="checkbox"/>		
	Diversity and richness (species count)	<input type="checkbox"/>		
	Microbial markers out of range - acetate, B. fragilis toxin, branched-chain amino acids (BCAA), beta-glucuronidase, butyrate, hexa-acylated lipopolysaccharide (hexa-LPS), hydrogen sulphide, 3-indolepropionic acid (IPA), methane, oxalate, propionate, trimethylamine (TMA)	<input type="checkbox"/>		
	Species table for more advanced users - to learn more visit Co-Education	<input type="checkbox"/>		
2 APPLY	Findings & insights	Based on patient symptoms, bowel habits, health history, allergies, intolerances, tolerability, goals, motivations	<input type="checkbox"/>	
		Prioritise insights based on the markers that need addressing the most (via results range or via health categories)	<input type="checkbox"/>	
		If there are no markers out of range, work on healthy microbiome foundations to help the patient improve their microbiome potential	<input type="checkbox"/>	
3 ADAPT	Treatment based on patient response & re-test results	Request further pathology or investigative testing, if necessary	<input type="checkbox"/>	
		Regular patient check-ins to monitor progress, compliance and treatment tolerability	<input type="checkbox"/>	
		Re-test between 3-6 months to assess treatment success	<input type="checkbox"/>	
		Maintain microbiome health	<input type="checkbox"/>	



6

# 1. ASSESS

## The Fundamentals of MetaXplore™ Report Interpretation

7

### MetaXplore™ Range Report Interpretation Checklist



1 ASSESS	Red flags (refer to a medical specialist, if necessary)	Faecal occult blood detected*	<input type="checkbox"/>	
		Calprotectin above 100 µg/g*	<input type="checkbox"/>	
		Lactoferrin above 7.2 µg/g*	<input type="checkbox"/>	
		Pancreatic elastase below 100 µg/mL*	<input type="checkbox"/>	
		Pathogens detected on diagnostic targeted pathogen panel**	<input type="checkbox"/>	
		Potential pathogens identified in metagenomic species table (search pathogen)	<input type="checkbox"/>	
	Gut terrain	Faecal pH*	<input type="checkbox"/>	
		Secretory IgA*	<input type="checkbox"/>	
		Zonulin*	<input type="checkbox"/>	
		Mucin degradation	<input type="checkbox"/>	
	Dysbiosis	Oral species	<input type="checkbox"/>	
		Diversity and richness (species count)	<input type="checkbox"/>	
		Microbial markers out of range - acetate, <i>B. fragilis</i> toxin, branched-chain amino acids (BCAA), beta-glucuronidase, butyrate, hexa-acylated lipopolysaccharide (hexa-LPS), hydrogen sulphide, 3-indolepropionic acid (IPA), methane, oxalate, propionate, trimethylamine (TMA)	<input type="checkbox"/>	
		Species table for more advanced users - to learn more visit Co-Education	<input type="checkbox"/>	
			<input type="checkbox"/>	



8

## Assess: Red flags

RED FLAGS

1. Faecal occult blood detected\*
2. Calprotectin above 100 µg/g\*
3. Lactoferrin above 7.2 µg/g\*
4. Pancreatic elastase below 100 µg/mL\*
5. Pathogens detected on diagnostic targeted pathogen panel\*\*
6. Potential pathogens identified in metagenomic species table

DETECTED

Faecal Occult Blood	
Calprotectin	<div style="display: flex; align-items: center;"> <div style="width: 100px; height: 10px; background: linear-gradient(to right, green, yellow, red);"></div> <div style="margin-left: 10px;"> <p>≤ 50.00 µg/g</p> <p style="color: red; font-weight: bold;">(110.99)</p> </div> </div>
Lactoferrin	<div style="display: flex; align-items: center;"> <div style="width: 100px; height: 10px; background: linear-gradient(to right, green, yellow, red);"></div> <div style="margin-left: 10px;"> <p>≤ 7.20 µg/g</p> <p style="color: red; font-weight: bold;">(12.87)</p> </div> </div>
Pancreatic Elastase	<div style="display: flex; align-items: center;"> <div style="width: 100px; height: 10px; background: linear-gradient(to right, red, yellow, green);"></div> <div style="margin-left: 10px;"> <p>≥ 200.00 µg/ml</p> <p style="color: red; font-weight: bold;">(62.13)</p> </div> </div>

DETECTED

● Proteobacteria
 Morganella morganii
0.51%

\*Available in MetaXplore GI & GI Plus only    \*\*Available in MetaXplore GI Plus only



9

## Assess: Why is it important to assess red flags first?

**Prioritise urgent health concerns:** Identifying red flags allows clinicians to address the most pressing health issues before focusing on broader patient goals, ensuring that severe conditions are managed first.

**Prevent misdiagnosis:** By assessing red flags early, clinicians can avoid overlooking significant health problems, reducing the risk of misdiagnosis or inappropriate patient management.

**Informed planning:** Understanding the presence of red flags informs the development of a more accurate and effective plans tailored to the patient's immediate needs.

**Referral to specialists:** Recognising red flags may indicate the need for referral to a specialist, ensuring that the patient receives the appropriate level of care for potentially serious conditions.




10

## Assess: When you need to refer on to a medical doctor

**Detection of serious pathogens:** If pathogens identified in the diagnostic panel or metagenomic species table indicate an **ACTIVE** infection with clinical symptoms, a referral is necessary for further investigation and management.

**Abnormal inflammatory markers:** Elevated markers such as calprotectin (>100 µg/g) or lactoferrin (>7.2 µg/g) suggest significant inflammation, which requires further evaluation.

**Digestive enzyme deficiency:** Pancreatic elastase levels below 50 µg/mL indicate possible exocrine pancreatic insufficiency, warranting a referral for specialised assessment.

**Presence of occult blood:** Detection of faecal occult blood necessitates immediate referral to rule out serious conditions such as gastrointestinal bleeding.


**Unexplained or severe symptoms:** If the patient's symptoms are severe, persistent or unexplained by the current findings, a referral is critical for comprehensive medical evaluation.



11


## Referral letter template

Available at [Co-Education](#)



12

## Assess: Gut terrain



\*Available in MetaXplore GI & GI Plus only \*\*Available in MetaXplore GI Plus only

CO-BIOME  
MICROBIA

13

## Assess: What are the signs of a compromised gut terrain?


**Altered faecal pH:** Changes in faecal pH can indicate an imbalance in the gut environment.

**Low/high secretory IgA levels:** Out-of-range secretory IgA suggests weakened mucosal immunity or activation, making the gut more vulnerable to infections and dysbiosis.

**Elevated zonulin levels:** High levels of zonulin are associated with increased intestinal permeability.

**Mucin degradation:** The presence of mucin-degrading species indicates compromised gut barrier function, which is essential for protecting the gut lining.

**Presence of oral species in the gut:** Finding oral species (>4) in the gut may suggest an altered gut environment.



CO-BIOME  
MICROBIA

14

## Assess: What does the gut terrain tell me about my patient?

### Gut transit time (faecal pH):

- Low pH might suggest rapid transit, which could lead to malabsorption or diarrhoea.
- High pH might indicate slow transit, possibly indicating constipation or increased fermentation in the gut. It can also help you understand SCFA absorption, crucial for maintaining gut health.

### Intestinal permeability (secretory IgA and zonulin):

- Elevated SIgA and zonulin levels may indicate a "leaky gut," which could lead to systemic immune activation or autoimmune conditions.
- Low SIgA might suggest a weakened mucosal defence, potentially leading to infections or poor response to pathogens.

### Mucosal protection (mucin):

- High levels of mucin-consuming microbes might signal an imbalance that could compromise the protective mucus barrier, leading to increased susceptibility to inflammation.

### Oral species:

- The presence of oral species in the gut may indicate translocation, suggesting compromised oral hygiene, stress, improper digestion or lack of stomach acid.



15

## Assess: Dysbiosis

### DYSBIOSIS

1. Diversity and richness (species count)

2. Microbial markers out of range –

- |  |  |
|--|--|
| <ul style="list-style-type: none"><li>• acetate</li><li>• <i>B. fragilis</i> toxin</li><li>• branched-chain amino acids (BCAA)</li><li>• beta-glucuronidase</li><li>• butyrate</li><li>• hexa-acylated lipopolysaccharide (hexa-LPS)</li></ul> | <ul style="list-style-type: none"><li>• hydrogen sulphide</li><li>• 3-indolepropionic acid (IPA)</li><li>• methane</li><li>• oxalate</li><li>• propionate</li><li>• trimethylamine (TMA)</li></ul> |
|--|--|

3. Species table



16




## Assess: What are the signs of functional dysbiosis?

**Diversity and richness imbalance:** Low species diversity and richness are signs of dysbiosis, potentially leading to decreased resilience of the gut microbiome and associated health issues.

**Imbalance of short-chain fatty acids (SCFAs):** Dysbiosis is often accompanied by an imbalance in SCFAs like butyrate, which is crucial for maintaining gut health and preventing inflammation.

**Overproduction of harmful compounds:** Overproduction of harmful compounds like hydrogen sulphide, hexa-acylated lipopolysaccharide (hexa-LPS), methane, and trimethylamine (TMA) indicates an imbalance in the gut environment, often linked to functional dysbiosis.

**Species overabundance:** An overabundance of certain microbial species, which can outcompete beneficial microbes, is a sign of functional dysbiosis.



CO-BIOME  
MICROBIA

17

## Assess: What does dysbiosis tell me about my patient?

**Microbial imbalance:** Dysbiosis indicates an imbalance in the gut microbiome, with a disruption in the balance of beneficial and harmful microbial species, which can impact overall health.

**Inflammatory potential:** The presence of dysbiosis suggests an increased potential for inflammation, both locally in the gut and systemically, potentially contributing to chronic conditions.

**Compromised gut barrier:** Dysbiosis may indicate a compromised gut barrier, increasing intestinal permeability and the risk of conditions such as IBS or autoimmune disorders.


**Need for targeted interventions:** Dysbiosis suggests the need for targeted dietary, lifestyle, or supplement interventions to restore balance and improve the patient's overall health.

CO-BIOME  
MICROBIA

18

## Whole microbiome approach

- It is important to consider both gut terrain and dysbiosis in the context of the whole microbiome, as well as your patient's symptom picture and case history.
- Which species are providing the functions you need? Rather than looking for one specific species, look at the functions being performed by the species present.
- Is your patient's diet and lifestyle impacting the potential for the microbiome to produce certain microbial markers?



CO-BIOME  
by MICROBA

19

## 2. APPLY

### Clinical Application of MetaXplore™ Report Interpretation

20

## MetaXplore™ Range Report Interpretation Checklist



<p>2 APPLY</p>	<p>Findings &amp; insights</p>	Based on patient symptoms, bowel habits, health history, allergies, intolerances, tolerability, goals, motivations	<input type="checkbox"/>	
		Prioritise insights based on the markers that need addressing the most (via results range or via health categories)	<input type="checkbox"/>	
		If there are no markers out of range, work on healthy microbiome foundations to help the patient improve their microbiome potential	<input type="checkbox"/>	
		Request further pathology or investigative testing, if necessary	<input type="checkbox"/>	



21

## How to prioritise insights for your patient management plan

### 1. Address red flags:

- Start by focusing on any red flags identified in the report. Provide insights or refer as required.

### 2. Prioritise managing gastrointestinal health markers that are out of range:

- Select insights that address gastrointestinal markers that fall outside the normal range.

### 3. Prioritise the microbial markers linked with the patient's clinical presentation:

- Address microbial markers where clinical symptoms are reported.

### 4. Identify any insights that can be applied across multiple markers:

- Look for insights that can positively impact several markers at once.


### 5. Consider which available insights the patient can follow/implement:

- Evaluate the practicality and feasibility of the insights from the patient's perspective. Consider factors such as the patient's lifestyle, preferences and ability to adhere to the recommendations.



22

## Support patient behaviour change with patient handouts and recipes



CO-BIOME  
-MICROBA

[Access Education via the Practitioner Portal](#)

23

## Common roadblocks to avoid when interpreting your MetaXplore™ report

**Getting bogged down in the species table:** Overfocusing on the species table can lead to missing the broader functional insights that are more clinically relevant.

**Looking at the microbiome as individual parts:** Treating the microbiome as separate, individual species rather than as an interconnected ecosystem can lead to a fragmented understanding.

**Focusing on just the celebrity species:** Placing too much emphasis on well-known or "celebrity" species instead of understanding the overall function of the microbiome.

**Searching for "missing" species:** Expecting specific species to be present in every microbiome and overlooking the fact that if a species isn't listed, it's not part of the patient's microbiome profile.

**Neglecting patient history:** Not integrating the patient's full medical history, including previous treatments and conditions, which can influence the interpretation of the report.

**Ignoring lifestyle factors:** Failing to account for lifestyle factors such as diet, stress and sleep that can significantly influence microbiome composition.

CO-BIOME  
-MICROBA

24

# Live Demonstration

**CO-BIOME**

Kylie Kingston  
Sex: Female  
DOB: 28 May 1968

**OVERVIEW**

- Summary
- Results
- Pathogen Panel

**HEALTH CATEGORIES**

- Intestinal Motility
- Intestinal Inflammation
- Intestinal Barrier
- Systemic Inflammation
- Detox / Retox
- Digestive Secretions

**OTHER TESTING**

- Emerging Metabolites
- Diversity

**DEEP ANALYSIS**

- Species Table
- Insights

**Targeted Pathogen Panel**

- Bacterial: NOT DETECTED
- Parasitic: NOT DETECTED

**Gastrointestinal Health Markers**

We've detected an out of range reading for 2 markers

- Calprotectin: IN RANGE
- Secretory IgA: OUT OF RANGE
- Pancreatic Elastase: BORDERLINE
- Faecal Occult Blood: IN RANGE
- Lactoferrin: IN RANGE
- Zonulin: OUT OF RANGE

**Health Categories**

- Intestinal Motility: 5 dots (3 green, 2 red)
- Intestinal Inflammation: 5 dots (3 green, 2 red)
- Intestinal Barrier: 5 dots (3 green, 2 red)
- Systemic Inflammation: 5 dots (all green)
- Detox / Retox: 2 dots (all green)
- Digestive Secretions: 1 dot (all orange)
- Emerging Metabolites: 5 dots (3 green, 2 orange)

**Microbiome Health**



25

**COMING SOON**  
MetaXplore™ Expert Summary

# Your fast track to microbiome interpretation

- A clear summary of results, right at the top.
- Understand indications quickly and save time.
- Informed by cutting-edge Microba science.

**CO-BIOME**

Kylie Kingston  
Sex: Female  
DOB: 28 May 1968

Report ID: R124080001261  
Sampling: 10:30am, 18 July 2023 (AEST) - Report released: 04:23pm, 28 August 2024 (AEST)

**Report Notes**  
Approved by Prof. Paul Griffin, Microbiologist.

**Expert Summary**

A clinical analysis of the patient's test results, summarised by our team of specialists.

There is evidence of intestinal permeability and an indication of reduced mucosal immune defence.

Zonulin is high (800µg/mL). Zonulin is a marker of increased intestinal permeability. Elevated zonulin is seen in patients with active coeliac disease, type 1 diabetes mellitus, metabolic syndrome, obesity, autoimmune disease, inflammatory diseases, neoplastic diseases, high faecal histamine, following high intensity exercise, and acute psychological stress.

Secretory IgA levels are low (244.83 µg/g). Secretory IgA plays a role in preventing adherence of microbes to mucosal sites. Reduced secretory IgA may be seen with increased fasting blood glucose in obese patients. Results should be interpreted in the context of the patient's clinical presentation.

Borderline low pancreatic elastase (187.19 µg/mL) should be considered in the context of clinical presentation. Review clinical symptoms and, if warranted, consider further investigation.

This microbiome contains methane-producing archaea at high levels compared to the healthy cohort. Around one-third of microbiomes contain methanogens. High levels of methane-producing archaea may contribute to reduced intestinal motility.

There is a low level of acetate producing microbes compared to the healthy cohort. Acetate will reduce inflammation when the system is in its normal state but will enhance the immune response when the immune system is activated.

The targeted pathogen panel did not detect the presence of pathogenic bacteria or parasites.

\*Some clinical microbiome measures are unvalidated and information provided is based on the professional expertise of the author.



26

1
ASSESS

Cutting down interpretation time to just minutes

### Expert Summary

A clinical analysis of the patient's test results, summarised by our team of specialists.

<p>There is evidence of intestinal permeability and an indication of reduced mucosal immune defence.</p> <p>Zonulin is high (800µg/mL). Zonulin is a marker of increased intestinal permeability. Elevated zonulin is seen in patients with active coeliac disease, type 1 diabetes mellitus, metabolic syndrome, obesity, autoimmune disease, inflammatory diseases, neoplastic diseases, high faecal histamine, following high intensity exercise, and acute psychological stress.</p> <p>Secretory IgA levels are low (244.83 µg/g). Secretory IgA plays a role in preventing adherence of microbes to mucosal sites. Reduced secretory IgA may be seen with increased fasting blood glucose in obese patients. Results should be interpreted in the context of the patient's clinical presentation.</p> <p>Borderline low pancreatic elastase (187.19 µg/mL) should be considered in the context of clinical presentation. Review clinical symptoms and, if warranted, consider further investigation.</p>	<span style="font-size: 24px;">✓</span> Gut terrain
<p>This microbiome contains methane-producing archaea at high levels compared to the healthy cohort. Around one-third of microbiomes contain methanogens. High levels of methane-producing archaea may contribute to reduced intestinal motility.</p> <p>There is a low level of acetate producing microbes compared to the healthy cohort. Acetate will reduce inflammation when the system is in its normal state but will enhance the immune response when the immune system is activated.</p>	<span style="font-size: 24px;">✓</span> Dysbiosis
<p>The targeted pathogen panel did not detect the presence of pathogenic bacteria or parasites.</p>	<span style="font-size: 24px;">✓</span> Red flags

\* Some clinical microbiome measures are unvalidated and information provided is based on the professional expertise of the author.

27

1  
**ASSESS**

➔

2  
**APPLY**

### Expert Summary

A clinical analysis of the patient's test results, summarised by our team of specialists.

There is evidence of intestinal permeability and an indication of reduced mucosal immune defence.

Zonulin is high (800µg/mL). Zonulin is a marker of increased intestinal permeability. Elevated zonulin is seen in patients with active coeliac disease, type 1 diabetes mellitus, metabolic syndrome, obesity, autoimmune disease, inflammatory diseases, neoplastic diseases, high faecal histamine, following high intensity exercise, and acute psychological stress.

Secretory IgA levels are low (244.83 µg/g). Secretory IgA plays a role in preventing adherence of microbes to mucosal sites. Reduced secretory IgA may be seen with increased fasting blood glucose in obese patients. Results should be interpreted in the context of the patient's clinical presentation.

**CO-BIOME**  
by MICROBA

### Insights

All Markers ▾
All Recommendation Types ▾
All Evidence Types ▾
All Evidence Grades ▾

✓IBS/food reactivity

CLINICAL INSIGHT

Offer glutamine supplementation to reduce intestinal permeability.

HUMAN B

SHOW REFERENCES

CLINICAL INSIGHT

Consider zinc carnosine supplementation to reduce intestinal permeability.

HUMAN C

SHOW REFERENCES

✓Constipation

RESEARCH INSIGHT

The probiotic *Lactobacillus reuteri* DSM 17938 may reduce methane production.

HUMAN D

SHOW REFERENCES

RESEARCH INSIGHT

Resistant starch supplementation has been shown to increase methane production. When aiming to reduce methane production, limiting or avoiding resistant starch supplementation may be effective.

HUMAN E

SHOW REFERENCES

✓Microbiome health

RESEARCH INSIGHT

In vitro stool fermentation studies suggest pectin may increase acetate production.

IN-VITRO PP, IV

SHOW REFERENCES

28

14

**“The new Expert Summary is a game changer, really insightful and can reduce the time taken on treatment decisions.”**



— *Practicing Naturopath*

*In a recent user testing session*

**Register to join the Early Access Program, today!**



## Key takeaways



## Live Q&A

Dr Brad Leech

Complete our survey and register for  
the MetaXplore™ Expert Summary  
Early Access Program



<https://t.maze.co/278981197>